

*Standard range*  
**ORTHOPEDIC  
PROSTHESES**



## Compliant with all ranges of knees and hips prosthesis

- Blisters for femoral knees, tibial plates
- Blisters for heads and cups
- Blisters for staples
- Blisters for femoral stem



**+** Safety

**+** No touch



### **+** Other products advantages

- Tray open : prosthesis displayed ready to use, this design improves the chirurgical acte.
- Tray closed : optimal protection for the prosthesis against impacts, falls or during the transport.
- Design adapted for heavy prosthesis : a raised bottom to absorb the kinetic energy during impact.

## Blisters for femoral hips

| Stem Revision | Reference | Stem position | Stem size (in mm) | Dimension outer blister (in mm) | Dimension inner blister (en mm) | Usefull dimension (in mm) |
|---------------|-----------|---------------|-------------------|---------------------------------|---------------------------------|---------------------------|
|               | K01       | 1             | <200              | 320 x 140 x 50                  | 298 x 118 x 43                  | 240 x 60 x 35             |
|               | K02       | 3             | <200              | 320 x 140 x 50                  | 298 x 118 x 43                  | 240 x 60 x 35             |
|               | K03       | 3             | <320              | 430 x 140 x 50                  | 408 x 118 x 40                  | 340 x 70 x 35             |

| Thomson and Moore prosthesis | Reference | Dimension outer blister (in mm) | Dimension inner blister (en mm) | Usefull dimension (in mm) |
|------------------------------|-----------|---------------------------------|---------------------------------|---------------------------|
|                              | K07       | 245 x 120 x 76                  | 223 x 98 x 67                   | 170 x 54 x 54             |
|                              | K12       | 290 x 140 x 81                  | 268 x 118 x 74                  | 230 x 62 x 62             |

## Blisters for heads and cups

| Reference | Dimension outer blister (in mm) | Dimension inner blister (en mm) | Usefull dimension (in mm) |
|-----------|---------------------------------|---------------------------------|---------------------------|
| K04       | 140 x 140 x 53                  | 118 x 118 x 49                  | 76 x 76 x 35              |
| K05       | 150 x 150 x 75                  | 128 x 128 x 73                  | 90 x 90 x 50              |

## Blisters for femoral knees and tibial plates

| Reference | Dimension outer blister (in mm) | Dimension inner blister (en mm) | Usefull dimension (in mm) |
|-----------|---------------------------------|---------------------------------|---------------------------|
| K06       | 150 x 150 x 76                  | 128 x 128 x 74                  | 75 x 75 x 55              |

## Blisters for staples and accessories

| Reference | Dimension outer blister (in mm) | Dimension inner blister (en mm) | Usefull dimension (in mm) |
|-----------|---------------------------------|---------------------------------|---------------------------|
| K08       | 98 x 88 x 54                    | 81,5 x 71,5 x 35                | 32 x 24 x 28              |
| K13       | 98 x 88 x 34                    | 81,5 x 71,5 x 25                | 32 x 24 x 23              |
| K1709     | 106 x 68 x 17                   | 88 x 50,5 x 15                  | 65 x 27 x 13              |
| K1709     | 106 x 68 x 17                   | Customized blister              |                           |



All products are packed in sealed double PE bags



Manufacturing in clean room ISO 7



Contact us

VELFOR GROUP

Route du Velay , 43500 SAINT-PAL-DE-CHALENCON, FR

Tel : +33 (0)4 71 61 30 51

Fax : +33 (0)4 71 61 38 44

contact@velfor.com

www.velfor.com

**velfor**  
Group  
Your development partner